

## MEDICAL SUPPLY & EQUIPMENT

Contact [eric@products4doctors.com](mailto:eric@products4doctors.com)  
cell 727-645-7774  
ofc. 727-999-3119



## MEDICAL SUPPLY & EQUIPMENT

Contact [eric@products4doctors.com](mailto:eric@products4doctors.com)  
cell 727-645-7774  
ofc. 727-999-3119

# ABOUT US

We treat our customers like family because we are a family owned and operated local business

## **We are a medical equipment superstore!**

We are one of largest retail supplier of medical equipment and supplies. with thousands of items in stock. We carry everything from compression stockings, braces, bath safety products, walkers, wheelchairs, CBD Gloves, needle & syringes, Diabetic or Blood Pressure Monitors, and so much more!

Our highly trained staff will help you make the appropriate choices for your needs, while providing you with friendly service and expert advice. Every effort is made to ensure that your experience is as pleasant and efficient as possible.



# DIRECTORY

#

RESPIRATORY CARE PRODUCTS

---

04

UROLOGY PRODUCTS

---

08

HYPODERMIC PRODUCT

---

19

MEDICAL EXAMINATION

---

26

WOUND DRESSING PRODUCTS

---

29

HOSPITAL UNIFORM

---

42



# RESPIRATORY CARE PRODUCTS



### YM-A001 SIMPLE OXYGEN MASK

- Made of medical grade PVC
- With adjustable nose clip
- With elastic strap
- With 2m tubing
- Size available S, M, L, XL
- Packed in individual PE Bag
- 100pcs/ctn.Dimension: 47\*39\*26cm

### YM-A002 VENTURI MASK

- Made of medical grade PVC
- With adjustable nose clip
- With elastic strap
- With 2m tubing
- Size available S, M, L, XL
- Packed in individual PE Bag
- 50pcs/ctn.Dimension: 45\*36\*20cm



### YM-A003 TRACHEOSTOMY MASK

- Made of medical grade PVC
- Flexible plastic mask with elastic head strap
- Excellent 360 angle rotation connector
- Excellent compatible with corrugated tubing
- Size available: Adult and Children
- Packed in individual PE Bag
- 100pcs/ctn.Dimension: 48\*36\*28cm

### YM-A002 VENTURI MASK

- Made of medical grade PVC
- With adjustable nose clip
- With elastic strap
- With 2m tubing
- Size available S, M, L, XL
- Packed in individual PE Bag
- 50pcs/ctn.Dimension: 45\*36\*20cm





### YM-A005 NEBULIZER WITH MOUTH MASK

- Made of medical grade PVC
- Used with any of the compressor/ nebulizers for treatment
- of asthmatic patients.
- Designed for comfortable administration of fine mists
- Of uniform density to reach deep into lung cavity.
- Simple three-piece design is easy to assemble and clean
- Break-resistant plastic provides the patient with 14 days of use
- With 15cm corrugated tube
- With 2m tubing
- With-piece and mouthpiece
- Packed in individual PE Bag
- 50pcs/ctn.Dimension: 49\*28\*22cm

### YM-A006 OXYGEN MASK WITH RESERVOIR BAG

- Made of medical grade PVC
- Reservoir bay with/without check valve
- One safety vent open.
- Mask with elastic strap
- Adjustable nose clip
- Non-rebreathing Mask provides a high concentration
- Of oxygen with 2m tubing.
- Sizes available: S, M, L, XL
- Packed in individual PE Bag
- 100pcs/ctn.Dimension: 47\*39\*25cm



### YM-A007 NASAL OXYGEN CANNULA

- Made of medical grade PVC
- Straight Tip, Flared Tip and Curved Tip is available.
- With 2m tubing
- Sizes available: S, M, L, XL
- Packed in individual PE Bag
- 100pcs/ctn. Dimension: 45\*38\*20cm





## YM-A007

### CONNECTING TUBE WITH YANKAUER HANDLE

- Made from non-toxic PVC, clear and soft
- Large lumen resists clogging and the transparency
- Allows clear visualization of liquids
- Crown tip, vented.
- With connecting tube
- Size: 1/4"\*1.8m, 1/4"\*3.6m, 3/16"\*1.8m, 3/16"\*3.6m,
- Packed in individual Blister Bag or Poly Bag
- 20pcs/ctn. Dimension: 38\*30\*25cm

## YM-A008

### ANESTHESIA BREATHING SYSTEM

- Adult/pediatrics; common/extension tube
- Construction with L connector, three-way connector.
- Breath air chamber, solenoid
- Function: connect the equipment, transfer the air to body
- Store claim: store in dark, dry and clean conditions,
- Do not give heavy press
- Form of packing: 30setjcarton
- Size of packing:
- Export carton: 660\*480\*375mm
- Gw: 5.8kgs



## YM-A009

### BACTERIAL/VIRAL FILTER

- Parts Shell membrane
- Material styrene
- Butadiene copolymer, Hydrophobic polypropylene
- 250pcs/60\*50\*50cm
- 8kgs/7.5kgs



# UROLOGY PRODUCTS



### YM-C001 2-WAY LATEX FOLEY CATHETER

- Made from natural latex
- Silicone coated
- Plastic valve or Rubber valve is available
- Color-coded for visualization of size Fr6-Fr 10( Pediatric)
- Length: 270mm Balloon Capacity: 3cc/Scc Fr12-Fr22(Standard)
- Length: 400mm Balloon Capacity: Scc/10cc/30cc
- Fr24-Fr26(Standard)
- Length: 400mm Balloon Capacity: 10cc/30cc
- Packed in individual Blister Bag
- 10Pcs/box, 500Pcs/ctn, Dimension: 52\*43\*47.Scm

### YM-C002 3-WAY LATEX FOLEY CATHETER

- Made from natural latex
- Silicone coated
- Plastic valve or Rubber valve is available
- Color-coded for visualization of size
- Fr16-Fr26(Standard)
- Length: 400mm Balloon Capacity: 30cc
- Packed in individual Blister Bag
- 10Pcs/box, 500Pcs/ctn, Dimension: 52\*43\*47.Scm




### YM-C003 2-WAY FEMALE FOLEY CATHETER

- Made from natural latex
- Silicone coated
- Plastic valve or Rubber valve is available
- Color-coded for visualization of size Fr12-Fr22
- Length: 240mm Balloon Capacity: 5cc/10cc/30cc Fr24-Fr26
- Length: 240mm Balloon Capacity: 10cc/30cc
- Packed in individual Blister Bag
- 10Pcs/box, 500Pcs/ctn, Dimension: 52\*43\*47.Scm



### YM-C004 2-WAY TIEMANN FOLEY CATHETER

- Made from natural latex
- Silicone coated
- Plastic valve or Rubber valve is available
- Color-coded for visualization of size Fr12-Fr22
- Length: 400mm Balloon Capacity: Scc/10cc/30cc Fr24
- Length: 400mm Balloon Capacity: 10cc/30cc
- Packed in individual Blister Bag
- 10Pcs/box, 500Pcs/ctn, Dimension: 52\*43\*47.5cm






### YM-C006 3-WAY DOUBLE BALLOONS FOLEY CATHETER

- Made from natural latex
- Silicone coated
- Color-coded for visualization of size Fr16-Fr24
- Length:400mm Balloon Capacity: 30cc(front)/50cc(rear)
- Packed in individual Blister Bag
- 10Pcs/box, 500Pcs/ctn, Dimension: 52\*43\*47.5cm

### YM-C007 2-WAY ALL SILICONE FOLEY CATHETER

- Made from high quality 100% silicone, medical grade
- With X-ray or without is available Fr8-Fr10(Pediatric)
- Length: 310mm Balloon Capacity: 3cc/5cc Fr12-Fr14(Standard)
- Length: 400mm Balloon Capacity:5cc/10cc Fr16-Fr24(Standard)
- Length: 400mm Balloon Capacity: 5cc/10cc/30cc
- Packed in individual Blister Bag
- 10Pcs/box, 500Pcs/ctn, Dimension: 52\*43\*47.5cm




### YM-C008 4-WAY DOUBLE BALLOONS FOLEY CATHETER

- Made from natural latex
- Silicone coated
- Color- coded for visualization of size Fr16-Fr24
- Length: 400mm Balloon Capacity: 30cc(front)/50cc(rear)
- Packed in individual Blister Bag
- 10Pcs/box, 500Pcs/ctn, Dimension: 52\*43\*47.5cm

### YM-C009 URETHRAL CATHETER(REL LATEX)

- Made from natural latex
- Silicone coated
- Fr6-Fr24
- Length: 400mm
- Packed in individual Blister Bag
- 10Pcs/box, 500Pcs/ctn, Dimension: 52\*43\*47.5cm





### YM-C010 MALECOT CATHETER

- Made from natural latex
- Silicone coated
- Fr14-Fr36
- Length: 400mm
- Packed in individual Blister Bag
- 10pcs/box, 500pcs/ctn, Dimension: 52\*43\*47.5cm

### YM-C011 ECONOMIC URINARY DRAINAGE BAG

- Capacity: 2000ml
- Without bottom outlet; with 90cm inlet tube
- With non-return valve
- Sterilization: EO
- Packing: 1pc/PE bag, 250pcstcarton
- Measurement of carton: 52\*32\*36cm



### YM-C012 ECONOMIC URINARY DRAINAGE BAG

- Capacity: 2000ml
- With bottom outlet; with 90cm inlet tube
- With non-return valve
- Sterilization: EO
- Packing: 1pc/PE bag, 250pcstcarton
- Measurement of carton: 540\*400\*270mm



### YM-C013 ECONOMIC URINARY DRAINAGE BAG

- Capacity: 2000ml
- With bottom outlet; with 90cm inlet tube
- With non-return valve, with T valve
- Sterilization: EO
- Packing: 1pc/PE bag, 250pcs/carton
- Measurement of carton: 52\*37\*31cm





### YM-C014 LUXURY URINARY DRAINAGE BAG

- Capacity: 2000ml
- Tape: square bottom, latex free, with pull push valve, with white sample port, with anti-reflux valve, with plastic tie
- Inlet tube: O.D 10mm; 100cm length
- Sterilization: EO
- Packing: 1pc/PE or blister pack, 100pcs/carton
- Measurement of carton: 54\*36\*32cm

### YM-C015 LUXURY URINARY DRAINAGE BAG

- Capacity: 2000ml
- Tape: square bottom, latex free With T valve; with white sample port With anti-reflux valve; without plastic tie
- Inlet tube: O.D 10mm; 100cm length
- Sterilization: EO
- From of Packing: 1pc/Biister pack, 100pcsjcarton
- Measurement of carton: 540\*360\*320mm



### YM-C016 LUXURY URINARY DRAINAGE BAG

- Capacity: 2000ml
- Tape: square bottom, latex free With T valve; with white sample port; With anti-reflux valve; with plastic tie
- Inlet tube: O.D 10mm; 100cm length
- Sterilization: EO
- From of Packing: 1pc/Biister pack, 100pcs/carton
- Measurement of carton: 540\*360\*320mm



### YM-C017 LUXURY URINARY DRAINAGE BAG

- Capacity: 2000ml
- Tape: square bottom, latex free with Single reversal valve; with needle sample port
- with anti-reflux valve with double plastic hanger and single rope with blue bed sheet clip
- Inlet tube: O.D 10mm; 100cm length
- Sterilization: EO
- From of Packing: 1pc/Biister pack, 40pcs/carton
- Measurement of carton: 430\*320\*380mm





### YM-C018 URINARY DRAINAGE BAG

- Capacity: 2000ml
- Tape: pyriform; latex free With T valve; with needle sample port; with anti-reflux valve; with double plastic hanger and rope with white bed sheet clamp
- Inlet tube: O.D 10mm; 110cm length
- Sterilization: EO
- From of Packing: 1pc/Biister pack, 40pcs/carton
- Measurement of carton: 430\*300\*380mm

### YM-C019 URINARY DRAINAGE BAG

- Capacity: 2000ml
- Tape: pyriform; latex free With T valve; with needle sample port; with anti-reflux valve with blue hanger and rope with white bed sheet clamp
- Inlet tube: O.D 10mm; 100cm length
- Sterilization: EO
- From of Packing: 1pc/Biister pack, 20pcs/carton
- Measurement of carton: 310\* 260\* 220mm



### YM-C020 URINARY DRAINAGE BAG

- Capacity: 2000ml
- Tape: pyriform; latex free With T valve; with needle sample port; with anti-reflux valve with white double plastic hanger and rope with white bed sheet clamp
- Inlet tube: O.D 10mm; 110cm length
- Sterilization: EO
- From of Packing: 1pc/ Biister pack, 40pcs/carton
- Measurement of carton: 430\* 300\* 380mm



### YM-C021 LUXURY URINARY DRAINAGE BAG

- Capacity: 4000ml
- with T valve; Latex free with anti-reflux drip chamber( one part) with needle sample port and tubing clip is optional with anti-reflux valve with reinforced double plastic hanger and rope with blue bed sheet clamp
- Inlet tube: O.D 11mm; 110cm length
- Sterilization: EO
- From of Packing: 1pc/Biister pack, 40pcs/carton
- Measurement of carton: 430\* 300\* 380mm






### YM-C022 LUXURY URINARY DRAINAGE BAG

- Capacity: 2000ml
- with T valve; Latex free with anti-reflux drip chamber(Three part) with needle sample port and tubing clip is optional with anti-reflux valve with reinforced double plastic hanger and rope with blue bed sheet clamp
- Inlet tube: O.D 11mm; 110cm length
- Sterilization: EO
- From of Packing: 1pc/ Biister pack, 40pcs/carton
- Measurement of carton: 430\* 300\* 380mm

### YM-C023 LUXURY URINARY DRAINAGE BAG

- Capacity: 2000ml
- with T valve; Latex free with anti-reflux drip chamber(Three part) with needle sample port and tubing clip is optional with anti-reflux valve with reinforced double plastic hanger and rope with blue bed sheet clamp
- Inlet tube: O.D 11mm; 110cm length
- Sterilization: EO
- From of Packing: 1pc/Biister pack, 40pcs/carton
- Measurement of carton: 430\*300\*380mm



### YM-C024 URINARY METER DRAINAGE BAG (DOUBLE CHAMBERS)

- Easy to read molded plastic meters are calibrated up to 400 ml urine meter
- with T valve on urine meter with twist valve on urine meter with needle sample port; with anti-reflux valve with reinforced double plastic hanger and rope with white bed sheet clamp
- Inlet tube: O.D 11mm; 130cm length
- Sterilization: EO
- From of Packing: 1pc/ Biister pack, 10Pcs/carton
- Measurement of carton: 480\*410\*250mm

### YM-C025 URINARY METER DRAINAGE BAG

- Capacity: 2600ml+500ml
- Latex free with T valve on urine meter with twist valve on urine meter with needle sample port; with anti-reflux valve with reinforced double plastic hanger and rope with white bed sheet clamp
- Inlet tube: O.D 11mm; 130cm length
- Sterilization: EO
- From of Packing: 1pc/Biister pack, 10Pcs/carton
- Measurement of carton: 480\*410\*250mm





### YM-C026 LEG BAG

- Capacity: 500ml and 1000ml
- Type: ellipse shape With green twist valve and green tip cover with anti-reflux valve; with double high elastic belts
- Sterilization: EO
- From of Packing: 1pc/PE pack, 48pcs/carton
- Measurement of carton: 300\*260\*150mm

### YM-C027 LEG BAG

- Capacity: 600ml
- Type: Square With blue pull-push valve and blue tip cover With anti-reflux valve
- With double high elastic belts With 120cm tube
- Sterilization: EO
- From of Packing: 1pc/PE pack, 48pcs/carton
- Measurement of carton: 300\*260\*150mm



### YM-C028 WOUND DRAINAGE SYSTEM

- Parts: Container, two to the connector, drainage pipe, connecting pipe, needle, non-return valve, and so on
- Main raw materials: PVC and/or silicone rubber
- According to drainage pipes and containers of different
- Materials used can be divided into PP PS SS three types
- According to the capacity of different containers
- can be divided into 200ml, 400ml, 500ml and 800ml.
- This product is used for abdominal, chest, breast and other parts of the fluid, pus and blood drainage.
- 200ml 20pcs/43\*32\*19cm
- 400ml 12pcs/43\*32\*19cm
- 500ml 24pcs/55\*32\*48cm
- 800ml 12pcs/45\*32\*48cm

### YM-C027 DISPOSABLE PEDIATRIC URINE COLLECTOR

- Size: 100ml, 200ml
- Component: pectic paper, PVC bag, sponge
- Sterilization: EO
- From of Packing: 1pc/PE pack, 5000pcs/carton
- Measurement of carton: 550\*540\*430mm





## YM-C030 WOUND DRAINAGE SYSTEM

- Parts: Container, two to the connector, drainage pipe connecting pipe, needle, non-return valve, and so on
- Main raw materials: PVC and/or silicone rubber
- According to drainage pipes and containers of different
- Materials used can be divided into, PP, PS, SS three types
- According to the capacity of different containers
- can be divided into 200ml, 400ml, 500ml and 800ml
- This product is used for abdominal, chest, breast and other parts of the fluid, pus and blood drainage

## YM-C031 ONE-PIECE CLOSED POUCH

- **PE foam adhesive**
- Highly resistant CO-EX film
- Foam basic plate is adherent, hypoallergenic
- And it is not easy to roll up
- Soft, comfortable
- It is easy to use
- It is for patients with solid faces and regular excretion



- 200ml 20pcs/43\*32\*19cm
- 400ml 12pcs/43\*32\*19cm
- 500ml 24pcs/55\*32\*48cm
- 800ml 12pcs/45\*32\*48cm



## YM-C032 ONE-PIECE CLOSED POUCH

- **Hydrocolloid adhesive**
- Highly resistant CO-EX film
- Hydrocolloid basic plate has strong viscosity
- hypoallergenic, comfortable and friendly to the skin Soft, comfortable It is easy to use
- It is for patients with solid faces and regular excretion

## YM-C033 ONE-PIECE MINI CLOSED POUCH

- **Hydrocolloid adhesive**
- Hydrocolloid basic plate has strong viscosity, soft, friendly and hypoallergenic to the skin
- Double-side non wove fabric lining is soft and opaque It is in miniature size, and easy to use
- Designed for one time use.
- Simply discard once full





YM-C034

**ONE-PIECE DRAINABLE POUCH WITH IRRIGATION PORT AND CLAMP CLOSURE**

- **Hydrocolloid adhesive**
- Highly resistant CO-EX film
- Hydrocolloid basic plate has strong viscosity, it can stick on the skin for 5 to 7 days
- It is soft, hypoallergenic
- The active carbon filter has the odor removal function
- The drainable design is easy to empty and clean the pouch
- Convenient for operation

YM-C035  
**ONE-PIECE DRAINABLE POUCH WITH TWIST-TIE CLOSURE**

- **Hydrocolloid adhesive**
- Highly resistant CO-EX film
- Hydrocolloid basic plate has strong viscosity, it can stick on the skin for 5 to 7 days It is soft, hypoallergenic
- The active carbon filter has the odor removal function
- The drainable design is easy to empty and clean the pouch
- Convenient for operation



YM-C036  
**ONE-PIECE DRAINABLE POUCH (CHILD)**

- Hydrocolloid basic plate has strong viscosity, soft, friendly and hypoallergenic to the skin
- There is nonwoven fabric or PE pouched lining on one side
- It is for pediatric stoma

YM-C037  
**TWO-PIECE CLOSURE POUCH**

- **Embedded**
- Soft and comfortable
- Embedded ring
- In miniature size and easy to use






YM-C038

## TWO-PIECE DRAINABLE POUCH WITH CLAMP CLOSURE

- **Embedded**
- Soft and comfortable
- Non-woven fabric or PE puncturing film liner
- The buckle ring is easy to attach the pouch with the skin barrier
- Active carbon filter



# HYPODERMIC PRODUCT



### YM-D001 DISPOSABLE INFUSION SET

- **Material: PVC**
- Small Chamber plastic needle+flow regulator+ rubber gum bubble
- With Latex Tube
- With Filter
- 500pcs 8x38x38cm

### YM-D002 DISPOSABLE INFUSION SET

- Small Chamber+stainless needle+flow regulator+rubber gum bubble
- With Latex Tube
- With Filter
- 500 pcs/58x38x38cm



### YM-D003 DISPOSABLE INFUSION SET

- Small Chamber+ plastic needle+flow regulator
- With Y site
- With Filter
- With luer lock adapter
- 500pcs/58x38x38cm



### YM-D004 DISPOSABLE INFUSION SET (DARK COLOR)

- **Material: PVC**
- Small Chamber+ plastic needle+flow regulator+rubber gum bubble
- With Latex Tube
- With Filter
- 500pcs/58x38x38cm





**YM-D005  
DISPOSABLE INFUSION SET WITH REGULATOR**

- Part: Infusion set+ IV regulator
- 500pcs/58x38x38cm

**YM-D006  
BURmE INFUSION SET**

- With Burette: 100ml, 150ml
- With Ysite
- With Filter
- Luer lock adapter with protector
- 100pcs/65x38x33.5cm



**YM-D007  
DISPOSABLE BLOOD TRANSFUSION SET**



- Designed for transfusion of blood or blood component
- Moulded latex flush ball injection membrane for intermittent medication
- Approximately 20 drops/ml
- Clear, transparent and soft tubing manufactured from medical grade PVC
- With 200 micron filter and 18 G hypodermic needle
- Options available
- With or without Airvent
- With or without Needle
- With or without "y" site( Latex or Latex-free)
- Luer lock or luer slip
- Packing; PE package or paper-plastic package
- 500pcs/65x38x33. Scm

**YM-D008  
DISPOSABLE BLOOD TRANSFUSION SET WITH REGULATOR**

- Part: Blood transfusion+IV regulator
- 500pcs/65x38x33. Scm





**YM-D009  
STOPCOCK THREE WAY**

- Parts: Body(PC), Plug(PE); Screw cap(ABS)
- Cold and heat tests, and the products will not have the phenomenon of crack
- Test pressure: 5.5kg.
- Products Thread: 6:100 screw thread standard
- This in line with international GB standards
- 1000pcs/58\*38\*38cm

**YM-D010  
SAFETY AUTO DESTRUCT SYRINGE**

- **Materials: PVC**
- Parts: Barrei+Piston+Piunger+Ciip+Needle (or without needle)
- Luer Lock or Luer Slip
- 2mi3000PCS/59\*44\*45cm      • 5ml2400 PCS /56\*44\*46.5cm
- 3m l 3000 PCS /59\*44\*45cm      • 10ml 1200 PCS /56\*44\*46cm



**YM-D011  
TUBERCLE BACILLUS SYRINGE**

- Size: 1ml
- Luer lock or luer slip
- 3600/16. 5x48. 5x38cm

**YM-D012  
DISPOSABLE SYRINGE**

- **Materaii:PVC**
- Parts: Barrei+Piston+Piunger+Needle(or without needle)
- Parts: Barrei+Piunger+Needle(or without needle)
- Luer lock or luer Slip/2 parts or 3 parts
- 1ml 3000PCS/59\*44\*45cm/3ml 3000 PCS/59\*14\*45cm
- 5ml 2400 PCS/56\*44\*46.5cm/10ml1200 PCS/56\*44\*46.5cm
- 20ml800 PCS/59\*44\*44cm/60ml400 PCS/59\*44\*14cm





## YM-D012 IRRIGATION SYRINGE

- Size: 50/60ml, 100ml, 150ml
- 50pcs/35x26x24cm

## YM-D013 COMESTIC SYRINGE

- PS material for the Barrel
- PP materil for the Plunger
- Size: 1 m l, Jml



## YM-D014 INSULIN SYRINGE

- Size: 100u
- Luer lock or luer slip
- 3600 pcs/16. 5x48. 5x3 Be m



## YM-D015 DISPOSABLE NEEDLE

- Size: 18G-31G
- High quality stainless steel
- Individually packing: blister packaging or in bulk.
- Inner packing: 100pcs/box
- Outer packing: corrugated carton





## YM-D016 SPINAL NEEDLE

- Spinal needle
- Needle materials: stainless steel
- Needle hub: transparent PP
- Quincke Tip: 18G-27G, 38-22Dmm
- Pencil point: 22G-27G, 90-22Dmm
- Color coded hub for size identification
- Packed in individual Blister Bag
- IDDDpcs/ctn, Dimension: 42\*27\*52cm

## YM-D017 EPIDURAL NEEDLE

- Size: 15G-18G, 80/90mm
- Needle materials: stainless steel
- Needle hub: transparent PP
- Color coded hub for size identification
- Packed in individual Blister Bag
- IODOpcs/ctn, Dimension: 42\*27\*52cm



## YM-D018 SCALP VEIN SET

- Butterfly shaped wings for better handing and attachment to skin
- siliconed cannula for smooth penetration
- Tube length: 15-30cm
- Gauge 18g-26g available
- Color coded hub for size identification
- Sterile by EO
- 1pc/poly bag 100pcs/box, 2000pcstcarton
- Carton size: 56\*25.5\*37.5cm



## YM-D019 DISPOSABLE EXTENSION SET

- Integrated Link Three ways stopcock facilitates the simultaneous administration of additional fluid from side port.
- Fitted with male luer lock at one end for safe and secure connection and the three way stop cock at the other end.
- Single Use, Sterile and non-Pyrogenic.
- Packed in individual Blister Bag

### • Tube Length:

- 10-25cm 50pcs/ctn, Dimension: 72\*56\*32cm
- 30-50cm 50pcs/ctn, Dimension: 72\*56\*32cm
- 8Dcm 50pcs/ctn, Dimension: 72\*56\*32cm
- 100-120cm 40pcs/ctn, Dimension: 72\*56\*32cm





## YM-D020 A.V.NEEDLE

- Mult size: 15G, 16G, 17G, 18G
- Needle with bevel edge with internal And external siliconisation, Ultra thin wall Stainless steel needle with/without backeye
- Self highly kink resistant tube
- Color code for size identification
- 300pcs/48\*44\*25.5cm

## YM-D021 I.V.CANNULA

- Needle size: 14,16,18,20,22,24G
- I.V.Cannula with injection port and suturable wings
- I.V.Cannula with suturable wings
- I.V.Cannula with injection port and without Wings
- Option Available:
- -PTFE/FEP/PU Flex Catheter
- -Hydrophobic filter
- -Clear or Radio Opaque Catheter
- Advantage of using PU Flex Catheter:
- -Kink free
- -The Catheter softens once inserted in the body
- -The properties of this Catheter are similar to PU (Polyurethane)



### FEATURES:

- Japanese Needle with 3-facet Needle point to reduce penetration force kink resistant and specially tapered
- Catheter for easy vein puncture with minimum trauma
- Transparent flash back chamber for easy visualization of blood to confirm venipuncture With luer lock plug
- Colour coding for each gauge as per international standards

## YM-D022 SAFEIY I.V.CANNULA

- Size: 14, 16, 18, 20, 22, 24G
- 1pc/blister bag, 1000pcs/ctn



## YM-D022 BLOOD BAG

Spedfication Beg	Capacity	Anticoagulant	Al foil pack/carton
Single	250ml-500ml	CPDA/CPO	6/120
Double	250mi-500ml	CPD,AJCPD	4/80
Triple	250mi-500ml	CPD,AJCPD+SAG-M	3/54
Quadruple	250mi-500ml	CPDA/CPD+SAG-M	3/45





# MEDICAL EXAMINATION



## YM-E001 SURGICAL GLOVES

- Description sterilized Surgical Gloves(Sterile By Gamma Ray).
- Powdered or powder-free
- Material: made of Natural Latex.
- Products Technical Datas:
- Specification (mm)

Size:	6	6.5	7	7.5	8	8.5	9
Length(min.mm)	260	260	270	270	270	280	280
Width(min,mm)	77±5	83±5	89±5	95±5	105±6	108±6	114±5
(70' cx168hours) tensile strength(n)	Before Again				After Again		
	>12				9		

Packing: SOPairs/box SOOPairs/ctn, 67\*26\*40cm, 15/13 Kgs

## YM-E002 LATEX EXAMINATION GLOVES

- Size: Large, Medium, Small
- Description: Powdered Or Powder-free
- Material: Made Of Natural Latex
- Products Technical Datas:
- Specification(mm):

Size:	XL	L	M	S	±10mm
Length(min.mm)	230				
Width(min,mm)	110	100	95	80	
Thickness Single (mm)	>0.08				
Physical Properties (mm)	26pcs sampled				
(100°C *22hours) tensile strength(mpa)	Before Again		After Again		
	21.0		16.0		
Ultimate elongation (%)	700		500		

Packing: IOOpcs/box, 2000pcs/53\*27\*32cm 18/16kgs





## YM-E003 NITRILE EXAMINATION GLOVES

- Description: nitrile Examination Gloves.
- Size: XL, Large, Medium, Small
- Description: Powdered or Powder free
- Products Technical Datas:
- Specification(m m):

Size:	XL	L	M	S	±10mm
Length(min.mm)	230				
Width(min,mm)	110	100	95	80	
Thickness Single (mm)	>0.08				
Physical Properties (mm)	26pcs sampled				
(100°C *22hours) tensile strength(mpa)	Before Again		After Again		
	14.0		13.0		
Ultimate elongation (%)	520		410		

Packing: 100pcs/box, 2000pcs/53\*27\*36cm, 18/16kgs

## YM-E004 VINYL GLOVES

- Material: Vinyl
- Product Name: Vinyl Gloves
- Grade: AQL4.0
- General Specification:
- Size: X-Large
- Weight: 6.0g±0.2
- Length: 240±10mm
- The Width of The Center of the Palm: 11.4±0.3cm
- Thickness: 0.05±0.01mm
- Content of Powderful: 0.1-0.2g
- Size: Middle
- Weight: S.Og±0.2
- Length: 240±10mm
- The Width Of The Centre OfThe Palm: 10.8cm±0.3
- Thickness: 0.05±0.01mm
- Content Of Powderful: 0.1-0.2g
- Size: Small
- Weight: 4.5g±0.2
- Length: 240±10mm
- The Width Of The Centre OfThe Palm: 10.6±0.3cm
- Thickness: 0.05±0.01mm
- Content Of Powderful: 0.1-0.2g
- We Can Make The Other Weight As Following:

S:	3.0g	3.5g	4.0g Etc
M:	3.5g	4.0g	4.5g Etc
L:	4.0g	4.5g	5.0g Etc





# WOUND DRESSING PRODUCTS



## YM-F001 GAUZE SWABS

- Material: 100% Cotton
- Cotton Yarn: 40's\*40's
- Mesh: 18\*11,18\*14,19\*9,19\*11,19\*15,20\*12,24\*20,26\*18,28\*24
- or Other Mesh Per Customer's Demand.
- Size: 2"\*2", 3"\*3", 4"\*4", 4"\*8"
- or Other Size Per Customer's Demand.
- Ply: 4Piy, 8Piy, 12Piy, 16Piy, etc.
- With or Without X-Ray
- Edge: Folded or Unfolded
- Color: White, Blue, Green
- Sterile or Non-Sterile
- Sterilization Method: EO or Gamma
- Sterile Packing: 1Pc/Two Side of Paper Pouch or One Side Paper+
- One Side Transparent PE Film Pouch, 2Pc/Pouch, 5Pc/Pouch or
- Other Packing Per Customer's Demand
- Non-Sterile Packing: 100Pcs/Paper Pack or 100 Pcs/Polybag


## YM-F002 GAUZE SWABS

- Material: 100% Cotton
- Cotton Yarn: 40's\*40's
- Color: White Bleached, Chlorine Free, Non-Toxic,
- Non Abrasive, Latex Free, Clean And Dust Free
- Mesh: 18\*11,18\*14, 19\*9, 19\*11,19\*15,20\*12, 24\*20,26\*18,28\*24
- or Other Mesh Per Customer's Demand.
- Size:36"\*50m, 36"\*100yards, 36"\*100m, etc.
- Ply: 1Piy, 2Piy, 4Piy etc.
- With or Without X-Ray
- Packing: 1roll/blue Kraft Paper Pack, Then Each
- Roll in Nylon Bag



## YM-F003 PREWASHED LAPAROTOMY SPONGE

- Material: 100% Cotton
- Cotton Yarn: 40's\*40's
- Color: White Bleached ,Chlorine Free ,Non-Toxic
- Non Abrasive ,Latex Free ,Clean And Dust Free
- With X-ray Detectable or Without X-ray Detectable
- Mesh:18\*11,18\*14, 19\*9, 19\*11,19\*15, 20\*12
- 24\*20,26\*18,28\*24
- or Other Mesh Per Customer's Demand.
- Size: 22.5cm\*22.5cm, 30cm\*30cm, 45cm\*45cm,etc.
- or Other Size Per Customer's Demand.
- Ply:4Piy, 8Piy, 12Piy etc.
- Sterile or Non-Sterile
- Sterilization Method: EO or Gamma
- Sterile Packing: SPes/Two Side of Paper Pouch or One
- Side Paper+ One Side Transparent PE Film Pouch,
- Non-Sterile Packing: 20Pcs/ Poly bag or in Bulk



## YM-F004

### PARAFFIN IMPREGNATED GAUZE PAD

- Description features: Sterile Paraffin Gauze Dressings.
- Are made of 100% cotton gauze swabs impregnated
- with white soft paraffin.
- Their multi-trulle structures are indicated as a primary
- wound contact dressing for manufactured from highabsorbent Gauze Swabs which have folded edges and
- thus prevent loose fibres and shredding.
- Size: 10cm\*10cm etc.
- Packing: each piece in foil pouch,
- 10pouch/inner box, 100pcs/ctn.
- Dim: 62\*28\*32cm, G.W: 8.0kgs, N.W: 7.0kgs

## YM-F005

### WOW GAUZE BANDAGE

- Material: 100% Cotton
- Cotton Yarn: 40's\*40's
- Color: White Bleached, Chlorine Free, Non-Toxic
- Non Abrasive, Latex Free, Clean And Dust Free
- Mesh: 18\*11,18\*14,19\*9,19\*11, 19\*15,20\*12
- 24\*20,26\*18,28\*24
- or Other Mesh Per Customer's Demand
- Width Size:1"/2"/3"/4"/5"/6"
- Length Size: 3yds/4yds/5yds/6yds/10yds
- or Other Size
- Sterile or Non-Sterile
- Sterilization Method: EO or Gamma
- Sterile Packing: Each Roll in Transparent PE Airproof
- Pouch or Other Packing per customer's Demand
- Non-Sterile Packing: 1 doz/Blue Kraft Paper Packing



## YM-F006

### NON-WOVEN SWABS

- Introduction: it is made from 30% terylene and 70% glue
- all cut edges are folded or unfolded and stiched evenly
- with a steady loop attached firmly to one corner. They are
- most useful for bleeding control or tissue protection in the
- operation rooms. it is normally with/ without X-ray
- detectable threads or tapes and can be dyed in various
- colors. It meets USP and European Phamacopoeia
- Gram of material: 30G/M2, 40G/M2, etc.
- Size: 2"\*2", 3"\*3", 4"\*4", 4"\*8"
- or Other Size Per Customer's Demand.
- Ply: 4Piy, 8Piy, 12Piy, 16Piy, etc.
- Color: White, Blue, Green
- Sterile or Non-Sterile
- Sterilization Method: EO or Gamma
- Sterile Packing: 1Pc{Two Side of Paper Pouch or One Side
- Paper+One Side Transparent PE Film Pouch, 2Pc/Pouch, 5Pc/
- Pouch or Other Packing Per Customer's Demand
- Non-Sterile Packing: 100Pcs/Paper Pack or 100 Pcs/Polybag



## YM-F007 ABSORBENT COTTON WOOL

- Introduction: Made from 100% natural cotton after cleaning
- purifying and bleaching, which is used for cleaning the wound
- It can be made into roll.
- Grade of Cotton: Natural or Super Absorbent
- Weight/Roll: 50Gram, 100Gram, 125Gram, 150Gram, 200Gram
- 400Gram, 454Gram, 500Gram, 1000Gram, etc.
- Packing: Each roll in Blue kraft paper packing, the details
- dimension of carton as following

Specification	Packing	Dimension of CTN(cm)
50g	1roll/Blue Kraft paper, 500rolls/carton	50*50*51
100g	1roll/Blue kraft paper, 200rolls/carton	61*31.5*61.5
200g	1roll/Blue Kraft paper, 1000rolls/carton	61*49*49
250g	1roll/Blue kraft paper, 1000rolls/carton	61*49*49
400g	1roll/Blue kraft paper, 36rolls/carton	49*31*49
454g	1roll/Blue kraft paper, 36rolls/carton	57*31*57
500g	1roll/Blue Kraft paper, 36rolls/carton	57*31*57
1000g	1roll/Blue kraft paper, 25rolls/carton	63*31*63

## YM-F008 ZIG ZAG COTTON WOOL

- Introduction: Made from 100% natural cotton after cleaning
- purifying and bleaching, which is used for cleaning the wound
- It can be made into piece.
- Grade of Cotton: Natural or Super Absorbent
- Weight/Roll: 25Gram, 50Gram, 100Gram, 200Gram, 250Gram
- 500Gram, 1000Gram, etc.
- Packing: sealed poly bag or drawstring bag, self-sealing bag



Specification	Packing	Dimension of CTN(cm)
25g	400bags/ctn	58*52*51
50g	200bags/ctn	64*58*46
100g	100bags/ctn	63*46*46
200g	50bags/ctn	63*46*48
250g	30bags/ctn	70*42*48
500g	20bags/ctn	46*40*47
1000g	10bags/ctn	68*50*43



## YM-F009 COTTON BALL

- Introduction: Made from 100% natural cotton after cleaning.
- purifying and bleaching, which is used for cleaning the wound.
- It can be made into ball.
- Grade of Cotton: Natural or Super Absorbent
- Weight/Roll: 0.3Gram, 0.5Gram, 1.0Gram, 2.0Gram, 3.0Gram, 4.0Gram, 5.0Gram, etc.
- Packing: sealed poly bag packing.
- Sterile or Non-Sterile
- Sterilization Method: EO or Gamma

Specification	Packing	Dimension of CTN(cm)
0.3G/Piece(sterile)	50piece/bag, 660bags/ctn	64*58*46cm
0.3G/Piece(sterile)	100piece/bag, 330bags/ctn	64*58*46cm
0.3G/Piece (Non-sterile)	1000piece/bag, 33bags/ctn	64*58*46cm
0.3G/Piece (Non-sterile)	2000piece/bag, 16bags/ctn	64*58*46cm
0.5G/Piece(sterile)	50piece/bag, 400bags/ctn	64*58*46cm
0.5G/Piece(sterile)	100piece/bag, 200bags/ctn	64*58*46cm

0.5G/Piece (Non-sterile)	1000piece/bag, 20bags/ctn	64*58*46cm
0.5G/Piece (Non-sterile)	2000piece/bag, 10bags/ctn	64*58*46cm
2G/Piece(sterile)	50piece/bag, 100bags/ctn	64*58*46cm
2G/Piece(sterile)	100piece/bag, 50bags/ctn	64*58*46cm
2G/Piece (Non-sterile)	1000piece/bag, 5bags/ctn	64*58*48cm
2G/Piece (Non-sterile)	2000piece/bag, 1bags/ctn	46*40*38cm

## YM-F010 GAUZE BALL

- Specification: 40's 18\*14mesh or 19\*11 mesh or 19\*15 mesh or 24\*20mesh or 26\*18mesh etc .
- Are available .
- Size: Dia.3cm or Dia. 4cm or Dia. 5cm are available
- Shape: Ball shape or peanut shape are available
- Color: pure white, green, blue are available



## YM-F011 COTTON SWAB(PLASTIC STICK)

- Specification: Single Tip For Medical Use Material of
- The Plastic Stick: PP, PS(can Break It off)
- Diameter of The Plastic Stick: 2.5mm-4.0mm
- Length of The Plastic Stick: 75mm-400mm
- Diameter ofThe Cotton Tip: 4-15mm
- Length of The Cotton Tip: 15-50mm
- 3"100pcs/bag 10bags/box, 20boxes/ctn, 47\*37\*23cm
- 6"100pcs{bag 10bags{box, 20boxes{ctn, 43\*37\*37cm
- 6"50pcs/bag 20bags/box, 8boxes{ctn, 63\*37\*49cm
- 8"50pcs/bag 20bags/box, 8boxes/ctn, 63\*37\*58cm



## YM-F012 EYE PADS, OVAL FORM

- Material: Cotton filled
- Cover: Gauze or Non-woven
- Type: Self-adhesive or without adhesive
- Size: 3cm\*6cm, 8cm\*6cm, etc.
- Sterile or Non-Sterile is available

### Packing:

- Size: 3cm\*6cm, 1pc/pouch, 50pouch/inner box
- 20box{ctn, Dim: 46\*42.5\*20cm
- Size: 8cm\*6cm, 1pc/pouch, 50pouch/inner box
- 20box{ctn, Dim: 57.5\*46\*20cm
- Size: 8cm\*6cm, 25pcs/pack, 40pcs/ctn
- Dim: 42\*26\*30cm

## YM-F013 COSMENC COTTON SWABS

- Introduction: Made from 100% Super Absorbent cotton
- and 100% natural paper fibre, the material is very soft
- Purifying with High temperature, which is used
- for cleaning the wound or dressing
- Shape: Round, or Square.
- Type: Plain, lattice and flower pattern
- Concave-convex dot, etc.
- Size: 5.5cm\*5.5cm, 6cm\*6cm, 7.5cm\*5cm
- Packing: 80pcs/Sealed plastic bag, 100pcs/Sealed plastic bag, 100pcs/inner box: ..
- or other package per customer's demand
- Normal Package:
- Size: 6.0cm{Round}, 100pcs{plastic bag
- 100bag/ctn, Dim: 62\*32\*38cm
- Size: 6.0cm\*6.0cm{Square), 100pcs/plastic bag
- 100bag/ctn, Dim: 62\*32\*38cm



## YM-F014 SPANDEX ELASTIC BANDAGE

- Material: 85% cotton With 15% spandex
- Crepe or Plain
- Color: natural, Flesh or Bleach
- Weight: 75g/m<sup>2</sup>, 80g/m<sup>2</sup>, 85g/m<sup>2</sup>, 90g/m<sup>2</sup>, etc.
- Width: 2", 3", 4", 5", 6"etc.
- Length: 4m, 4.5m, 5m, 6m, etc.
- Sterile or Non-Sterile
- 1roll{poly Bag.12roll{bag}



### YM-F015 PBT ELASTIC BANDAGE

- Material: 65% cotton With 35% polyester
- Color: Bleach
- Weight: 28g/m', 29g/m', 30g/m' etc.
- Width: 2", 3", 4", 5", 6"etc.
- Length: 4m, 4.5m, 5m, 6m, etc.
- Sterile or Non-Sterile
- trol/poly Bag.12roll/bag

### YM-F016 COTTON ELASTIC BANDAGE

- Material: 100% cotton
- Color: Bleach
- Weight: 75g/m2, 80g/m2, 85g/m2, 90g/m2, etc.
- Width: 2", 3", 4", 5", 6"etc.
- Length: 4m, 4.5m, 5m, 6m, etc.
- Sterile or Non-Sterile
- 1roll{poly Bag.12roll/bag



### YM-F017 HIGH ELASTIC BANDAGE

- Material: 62% cotton With 38% polyester
- Color: Bleach, Natural
- Weight: 80g/m', 100g/m', etc.
- Width: 2", 3", 4", 5", 6"etc.
- Length: 4m, 4.5m, 5m, 6m, etc.
- Sterile or Non-Sterile
- 1roll/poly Bag.12roll/bag



### YM-F018 NET BANDAGE

- Material: Acrylic with Elastic
- Color: White
- Weight: 120g/m2, etc.
- **Size:** 0.8cm\*25m, 1.7cm\*25m, 2cm\*25m, 2.2cm\*25m, 2.6cm\*25m, 3.8cm\*25m, 4.4cm\*25m, 5.5cm\*25m, 6.8cm\*25m, 8cm\*25m





### YM-F019 NON-WOVEN SELF-ADHESIVE BANDAGE

- Material: Non-woven PP
- Color: White, blue green, red, brown, yellow etc.
- Width: 2", 3", 4", 5", 6" etc.
- Length: 4m, 4.5m, 5m, 6m, etc.
- Sterile or Non-Sterile
- 1roll/poly Bag.12roll/box

### YM-F020 COTTON SELF-ADHESIVE BANDAGE

- Material: Cotton
- Color: White ,blue green, red, brown, yellow etc.
- Weight: 75g/m<sup>2</sup>, 80g/m<sup>2</sup>, 85g/m<sup>2</sup>, etc.
- Width: 2", 3", 4", 5", 6", etc.
- Length: 4m, 4.5m, 5m, 6m, Etc.
- Sterile or Non-Sterile
- 1roll/poly Bag.12roll/box



### YM-F021 PLASTER OF PARIS BANDAGE

- Material: Cotton With Plaster
- Color: White
- Weight: 390g/m<sup>2</sup>, 400g/m<sup>2</sup>, 410g/m<sup>2</sup>, etc.
- Width: 2", 3", 4", 5", 6", etc.
- Length: 4m, 4.5m, 5m, 6m, etc.
- Sterile or Non-Sterile
- 1roll/poly Bag.12roll/box

### YM-F022 SELF-ADHESIVE ELASTIC BANDAGE

Material: Cotton With Spandex  
Crepe or Plain  
Color: White, Skin  
Weight: 75g/m<sup>2</sup>, 80g/m<sup>2</sup>, 85g/m<sup>2</sup>, 90g/m<sup>2</sup>, etc.  
Width: 2", 3", 4", 5", 6", etc.  
Length: 4m, 4.5m, 5m, 6m, etc.  
Sterile or Non-Sterile  
1roll/poly Bag.12roll/box





### YM-F023 FIRST- AID BANDAGE

- Cotton filled Size: L, M, S
- Packing: Vacuumize
- Packing Size: Small, 1pc/Vacuumize pack
- 400pack/ctn, Dim: 50\*44\*36cm
- Size: Large, 1pc/Vacuumize pack
- 400pack/ctn, Dim: 50\*46\*43.5cm

### YM-F024 TRIANGLE BANDAGE

- Material: 100% cotton or Rayon Fabric
- Cotton Yarn: 40\*34, 50\*30, 48\*48, etc.
- Size: 36"\*36"\*51", 40"\*40"\*56",
- With or without safety pin Colored: Grey or Bleached.
- Packing: 1pcfsealed poly bag. 240bagjctn.
- Dim: 45\*38\*22cm.
- G.W: 6.5kgs, N.W: 5.0kgs



### YM-F025 SPORTS TAPE

- Ideal for sport protection, water proof, permeability
- Different colors available
- Size available: 2.5cm\*10y, 3.8cm\*10y, 5cm\*10y
- 2.5cm\*15y, 3.8m\*15y, 5cm\*15y
- Packing 4roll/box, 300box/ctn
- Dimension: 34.5\*26.5\*28cm



### YM-F026 FIBERGLASS CASTING TAPE

- Material: Fiberglass  
Size: 5cm\*360cm, 7.5cm\*360cm,  
10cm\*360cm  
12.5cm\*360cm, 15cm\*360cm  
Packing: 1pc/Aiuminum Foil Bag  
50pcs/box  
100box/carton





## YM-F027 ICE BANDAGE

- Ice wrap is a cold bandage which uses a revolutionarily new technique for treating injuries, muscle pain and strains.
- This bandage reduces recovery time from injuries by combing both the cooling and compressing process in an easy-to-use way.
- A specially formulated solution is infused into the elastic bandage that begins to evaporate upon contacting skin, thereby causes cooling of the applied area. Lcewrap reduces swelling and provides sufficient cooling to make it a very effective treatment.

### Size available:

10cm*3.2m	120rolls/ctn	66*33*26cm
7.5cm*3.2m	120rolls/ctn	66*33*26cm

10cm*3.2m	48rolls/ctn	36*27*42cm
7.5cm*3.2m	80rolls/ctn	36*27*42cm


## YM-F028 ZINC OXIDE PLASTER

- Material: 10% Zinc oxide with 90% hot melt adhesive
- Color: White, Flesh.
- Weight of adhesive: 80-SSg/m'
- Width: 1.25cm, 2.5cm, 5cm, 7.5cm, 10cm, etc.
- Length: 5m, 10m, etc.
- Packing: Plastic can, Tinplate or simple packing



## YM-F029 NON-WOVEN SURGICAL TAPE

- Material: Non-woven with adhesive
- Color: White, Flesh
- Weight of adhesive: 30-32g/m'
- Width: 1.25cm, 2.5cm, 5cm, 7.5cm, 10cm, etc.
- Length: 5m, 10m, etc.
- Packing: Plastic can, Plastic cutter or simple packing



## YM-F030 PE SURGICAL TAPE

- Material: Polyethylene with adhesive
- Color: Transparent
- Weight of adhesive: 30-34g/m'
- Width: 1.25cm, 2.5cm, 5cm, 7.5cm, 10cm, etc.
- Length: 5m, 10m, etc.
- Packing: Plastic can, Plastic cutter or simple packing.



### YM-F031 SILK SURGICAL TAPE

- Material: Silk with adhesive
- Drilled or non-drilled is available
- Line shape or zig shape is available
- Color: White
- Weight of adhesive: drilled 50-54g/m'
- Non-drilled 32-34g/m'
- Width: 1.2Scm, 2.Scm, Scm, 7.5cm, 10cm, etc.
- Length: Sm, Sy, 10y, etc.
- Packing: Plastic can, simple packing


### YM-F032 WOUND ADHESIVE PLASTER

- Material: 100% cotton or Rayon Fabric
- Cotton Yarn: 40\*34, 50\*30, 48\*48, etc.
- Size: 36"\*36"\*51", 40"\*40"\*56",
- With or without safety pin Colored: Grey or Bleached.
- Packing: 1pcfsealed poly bag. 240bagjctn.
- Dim: 45\*38\*22cm.
- G.W: 6.5kgs, N.W: 5.0kgs




### YM-F033 DRILLED PLASTER

- Material: Drilled Zinc oxide with hot melt adhesive.
- Waterproof or non-water proof is available
- With drug or without is available
- Color: White, Flesh
- Weight of adhesive: 60-65g/m'
- Width: 2.5cm, scm, 7.5cm, 10cm, 12cm, etc.
- Length: 1m, Sm, etc.
- Packing 1roll/box, 80roll/carton



### YM-F034 HEALING CAPSICUM PLASTER

Capsicum plaster uses the full cotton cloth as the raw material, which was coated on capsicum consisted zinc oxide adhesive. Capsicum plaster can make the partial skin generate physical reaction to accelerate blood fluxion This plaster has drilled type and not-drilled type. All the products adopt the anti-sticky paper and film as the interleaver Can be used by one piece or cut size conveniently and easily.





### YM-F035 SELF-ADHESIVE WOUND DRESSING

- Material: PU, Non-woven, etc.
- Color: White, Flesh, etc.
- Weight of adhesive: 70-75g/m<sup>2</sup>
- Shape: Round corner, Square corner.
- With mat or without is available.
- Size: 5cm\*7cm, 6cm\*8cm, 10cm\*10cm, 10cm\*12cm, 10cm\*15cm, 10cm\*20cm, 10cm\*25cm, 10cm\*30cm, etc.
- Packing: 50pcs/box

### YM-F036 INFUSION PLASTER

- Material: Non-woven
- Color: White, Flesh
- Weight of adhesive: 30-32g/m<sup>2</sup>
- Shape: Rectangle
- Size: 75\*25cm, 75\*35cm, etc.
- Packing: 10Dpcs/box, 100boxcarton



### YM-F037 SELF-ADHESIVE EYE PAD

- Material: Non-woven
- Color: White, Flesh
- Weight of adhesive: 30-32g/m<sup>2</sup>
- Shape: Oval
- Size: 80cm\*50cm, 85cm\*55cm, 95cm\*65cm, etc.
- Packing: 50pcs/box, 100boxcarton



### YM-F038 IV CANNULA DRESSING

Material: Transparent Film Dressing  
Hydrocolloid I.V dressing  
PE Foam I.V dressing  
Non-woven I.V dressing





## YM-F043 ALCOHOL PREP PAD

- Contents: 70% Isopropyl Alcohol Non-woven
- Size: 65\*30mm, 65\*58mm, 140\*180mm or Other
- Size Per Customer's Demand
- Type: 1Piy, 2Piy, 4Piy, etc.
- Packing Size: 50\*50mm, 60\*85mm, etc.
- Packing Quality: 1Pc/Aluminium foil bag
- 200pcs/box, 20000pcs/ctn, 55\*49\*31cm

## YM-F044 POVIDONE-IODINE PREP PAD

- Contents: 10% Povidone-iodine Solution Nonwoven
- Size: 65\*30mm, 65\*58mm, or Other Size Per Customer's Demand.
- Type: 1Piy, 2Piy, 4Piy, etc.
- Packing Size: 50\*50mm, etc.
- Packing Quality: 1pc/Airtight Pack
- 200pcs/box, 20000pcs/ctn, 55\*49\*31cm.



## YM-F045 ADHESIVE TAPE REMOVER PAD

- Contents: Adhesive Tape Remover Pad
- Non-Woven Size: 65\*30mm
- Packing Size: 50\*50mm
- Packing Quality: 100Pcs/Box, 200Pcs/Box.
- 50boxes/carton, 53\*24\*29cm.



## YM-F046 STING RELIEF

Sting relief used in wiping the skin area where was stung by mosquito, bee and other insects for relief of pain and itch.

30\*60mm\*2non-woven pad wrapped in a foil paper sachet  
sized 50\*50mm

Dimension: 54\*23\*28cm





# PERSONAL PROTECTION PRODUCTS



## YM-G001 ROUND NURSE CAP

### Description:

- Manufactured From High Quality Non-woven Polypropylene Lightweight, Cool, Comfortable Material-Soft Elastic Headband Suited To Lab, Industry, Hospital And Daily Home Use

### Specification:

- Non Woven PP Material: 10g,12g,14g, 16g,18g
- Available in 18", 21" or 24" Sizes
- White, Green or Blue, Pink, Yellow, Red, etc.
- Sterile Method: None or EO Sterile
- Packing: 100pcs/Poly Bag 10Bags/Carton
- Carton size: 240\* 240\*340mm

## YM-G002 CLIP CAP

### Description:

- Manufactured From High Quality Non-woven Polypropylene High Quality Material For Increased Durability Convenient Dispenser Light weight, Cool, Comfortable Material-Soft Elastic Headband, and Comfortable Fold Type Suited To Lab Industry. Hospital And Daily Home Use

### Specification:

- Non Woven PP Material:
  - Sterile Method: None or EO Sterile
- 10g,12g, 14g, 16g, etc.
- Color: White, Green Blue, Pink, Red, Yellow
  - Packing: 100pcs/Poly Bag 10Bags/Carton
- Size: 20",21" or 24"
  - Carton size: 530\*240\*350mm
- Type: Single Elastic/Double Elastic
  - Carton size: 240\* 240\*340mm



## YM-G003 DOCTOR CAP

### Description:

- Manufactured From High Quality Non-woven Polypropylene High Quality Material For Increased Durability Convenient Dispenser Light weight, Cool, Comfortable Material-Soft Elastic Headband, and Comfortable Fold Type Suited To Lab Industry. Hospital And Daily Home Use

### Specification:

- Non Woven PP Material: 20g, 25g, 30g, 40g, etc.
- Color: Green/Blue
- Type: Tie or Elastic back
- Size: 61 \*14cm; 64\*23cm etc
- Latex free
- The SMS, Spun Lace Fabric is available
- Packing: 100pcs/Poly Bag 10Bags/Carton
- Size of packing:
- Export carton: 50\*50\* 30cm or 18\*30\*30cm
- Gross Weight: 7.0kgs
- Net Weight: 6.0kgs



## YM-G004 DOCTOR CAP

### Description:

- Manufactured From High Quality Non-woven Polypropylene High Quality Material For Increased Durability Convenient Dispenser Light weight, Cool, Comfortable Material-Soft Elastic Headband, and Comfortable Fold Type Suited To lab, Industry Hospital And Daily Home Use

### Specification:

- Nylon Net Fabric
- Color: Blue
- Type: Tie or Elastic back
- Size: 18",21"or 24"
- Packing: 100pcs/Poly Bag 10Bags/Carton
- Size of packing:
- Inner box: 280\*260\*100mm
- Export carton: 510\*290\*530mm
- Gross Weight: 6.0kgs
- Net Weight: 5.0kgs

## YM-G005 NON WOVEN FACE MASK

### Description:

- Conforms To American Nelson N95 Test Standards Mini-Fiber High Electrostatic Filter Filtration Efficiency of Dust and Powder is over95% Repel Continually Contamination of Dangerous Blood And Saliva Cover is Water Resistant. Polypropylene Cloth Low Breathing Resistance and High Filtration Efficiency

### Specification:

- Ear-Loop Type: 1ply, 2-ply or 3ply
- Detailed Weight: 1Piy-38G/m<sup>2</sup>, 2Piy: 25g/m<sup>2</sup>+18g/m<sup>2</sup>, 3Piy-25G/m<sup>2</sup>+20g/m<sup>2</sup> +16g/m<sup>2</sup>
- Size: 18\*9cm
- Color: White, Green, Blue, etc.
- Packing: 50pcs/Box, 40boxes/ Carton
- Carton size: 420\*390\*370mm



## YM-G006 NURSE PAPER FACE MASK

### Description:

- Our Most Economical Face Masks, High Quality WetStrength paper Comfortable Elastic Ear Loops for Easy Donning Product

### Specification:

- 2Piy or 1Piy
- Size: 20\*7cm
- Color: White
- Suggested Applications : Polluted Environment
- Especially Good for Absorbing and Filtrating The Exhausted Gas
- From Motor Vehicles.
- Packing: 100pcs/Box, 100boxes/Carton



## YM-G007 DUST-PROOF FACE MASK

### Specification:

- Conforms to FFPI and FFP2
- Fiberglass Free, odorless, Non-irritating
- To the Skin Adjustable Contoured Aluminum
- Nosepiece Helps To Ensure Proper Fit.
- Low Cost And Easy To Use
- Suggested Applications: Factory, Garden
- Process, Maintain.
- House Cleaning, Motorcycle Riders, Wood and Any Toxic Free Protection Non-woven PP
- Material: 16g, 18g, 20g, etc.
- Color: White, Blue, Dark Blue
- Packing: 50Pcs/inner Box, 20Boxes/Carton
- Size of packing:
- Inner box: 115\*115\*125mm
- Export carton: 600\*260\*250mm
- Gross Weight: 4.2kgs
- Net Weight: 3.6kgs

## YM-G008 FACE GUARD

### Specification:

- Conforms To American Nelson N95 Test Standards
- Ear-Loop Type: 3-Ply or 2-Ply Non Woven, Polypropylene Cloth, Low Breathing Resistance and High Filtration Efficiency. and PVC Transparent Shield Mini-Fiber High Electrostatic Filter, Filtration Efficiency of Dust and Powder is over 95%.
- Repel Continually Contamination of Dangerous Blood
- And Saliva Continually. Inner Cover is Well Sweet Absorbed, While Outer Cover Is water Resistant.
- Packing: 25Pcs/inner Box, 20boxes/Carton
- Size of packing:
- Export carton: 540\*340\*210mm
- Gross Weight: 4kgs
- Net Weight: 3kgs



## YM-G009 SURGICAL HOOD

### Description

- A Disposable Surgical Hood Which Conforms To The Wearer's Head .and is Usable Either with a traditional Surgical Mask Alone or With Such a Mask In Combination With a Laminar Air Flow Face Plate Contains Bacterial Fallout From The Wearer's Head and Helps Maintain Sterile Conditions in Hospital Operating Rooms.

### Specification:

- Non-Woven PP Material: 16g, 20g, 25g, 30g, etc.
- Size: 30\*24cm
- Color: White, Green, Blue
- Packing: 100pcs/poly bag; 10bags/Carton
- Size of packing:
- Export carton: 380\*300\*300mm
- Gross Weight: 7.5kgs
- Net Weight: 6.2kgs



## YM-G010 NEW LYPE DUST PROOF FACE MASK

### Description:

- Conforms To American Nelson N95 Test Standards Ear-Loop For Easy Use, 160g/180g/200g Low Breathing Resistance With High Filtration Efficiency Adjustable Nose Seal, Easy To Fold and Carry Special 3-Piy Design Respirator Produced From High Filtrable And Electrostatic Material With High Efficient Activated Carbon Cloth For Relief From The Irritating Effects of Nuisance Odor, Acid Gas, Hydrogen Fluoride and Ozone

### Suggested Applications:

- Polluted Environment,
- Especially Good For Absorbing and Filtrating The
- Exhausted Gas
- From Motor Vehicles.
- Packing: 500pcs/inner box; 20boxes/ Carton
- Size of packing:
- inner box: 115\*115\*128mm
- Size of outer box: 600\*270\*250mm
- Gross Weight: 5kgs Net Weight: 4kgs

## YM-G011 N95 FACE MASK(WITH/WITHOUT VALVE)

### Specification:

- Conforms to International Standards EN149: 2001 FFPI FFP2 CE0194 AS/NZS 1716: 1994 NIOSH N95 Low Breathing Resistance with High Filtration Efficiency Adopted The Polypropylene Material Which is Without Any Toxic Chemicals or Stimulant Adjustable Contoured Aluminum Nosepiece and a And Increase Worker Comfort. Easy-Breathing Exhalation Valve With Inner Silicon-Piece Can Facilitate The Air Exchange and increase Worker Comfort Special 3-Piy Design Respirator Produced From High Filtrable And Electrostatic Material For Resisting Normal Toxic and Tiny Micro Particle.

### Suggested Applications:

- Chemurgy, Spray-Paint/Printing,
- Dyeing, Garage, Factories of Electronic Aid and food
- Disinfecting and The Insalubrious Workshops
- From of Packing: 20pcs/inner box; 20boxes/Carton
- Size of packing:
- inner box: 120\*120\*260mm
- Export carton: 620\*260\*540mm
- Gross Weight: 8kgs Net Weight: 7kgs



## YM-G012 FILTERING HALF MASK

- Flat-folded CE EN 149: 2001
- Protection against solid and fluid particles
- (FFPI and FFP2)
- And protection against toxic particles(FFP3 only), single
- packed without valve, FFPI/with valve FFP2/with valve,
- FFP3
- Package: 20pcs/box; 400pcs/Carton
- Carton size: 57\*43\*32cm



## YM-G013 SURGEON GOWN

### Specification:

- Materials Non-Woven PP/PP+PE+PE/SMS
- Color White/blue/Yellow/Green
- Weight: 40/45/50GSM PP or SMS
- Size: S: 110\*154cm/ M: 120\*150cm /L: 130\*155cm/XL: 140\*160cm
- Suited to Hospital Use. 100pcs/Bag 100pcs/Carton With Ties At Back of The Neck And Waist
- Comfortable Sleeve Tailoring For Maximum Freedom of Movement
- Soft Kint Cuffs For Wrist Comfort Secure Sleeve Seams For Protection
- Extra Protection Reinforcement on sleeve and Chest and Full-body on selected Gowns
- Cloth Like Fabric Breathable Quiet and Supply With Knitted Cuff Vecro Neck Back Closure,
- Side Waist Band Reinforced Arm And Additional Liner at
- Front is Available
- !Hand Towel is Available
- Sterile Method :Non Sterile or EO Sterile
- Packing: 100pcs/poly bag or 1pc/ Sterile Bag
- SMS Surgical Gown
- Natural Latex Free, Water Repellent, Protection Against Splatter of Blood. Breathable, Soft
- 100% PP Spunbond+Meltblown(SMS)
- Anti-static is available
- Spunlace Surgical Gown
- Anti-Bacterial Alcohol-Repellent, Blood- Repellent Waterproof Dustproof Breathable, Soft Comfortable Mobility Made from Spun lace Non-Woven Packing: 1pc/Sterial Pack 25pcs/ Carton

## YM-G014 DISPOSABLE UNIFORM ISOLATION GOWN

### Specification:

- Materials: offer Economical
- Comfortable Reliable and Protection
- Non-Woven PP/PP+PE
- Full Length For Maximum Coverage, Latex Free
- Sewn Seams Provide Super Strength
- And Extra Length On Waist Ties Allow Gown
- To Be Easily Secured In Front
- With Ties at Back of The Neck and Waist
- Color: White/blue/Yellow/Green
- Weight: 18/25/30/40GSM or 18/25/30/40 GSM PP+15PE
- Size: 115\*137cm or 120\*140cm etc.
- Elastic Wrist or With Knitted Cuff
- 100pcs/bag, 100pcs/Carton
- OR 1pc/Sterile Pack; 50Pcs / Carton
- Sterile Method: EO Sterile,
- Suited To Lab, Industry, Hospital and Daily Home Use.
- Form or Packing:
- 100Pcs/Poly Bag or 1pc/Ster or Packing
- Packing: 1Pcs/Poly Bag, 10Bags/ Carton
- Carton size: 560\*320\*350mm





## YM-G015 CPE GOWN

### Specification:

- This kind of CPE gown can be used in hospital, laboratory Food industry, any where need protection

### Description:

- Prevent and isolate dust, particles, Alcohol, blood, bacteria and virus invading. Can be used as a visit coat Prevent te cross infection in medical treatment at sanitation field Economical and disposable.
- Accord with the hygienic standard

- Materials: CPE
- Characteristics: prevent and isolate dust, particle, alcohol, blood bacteria and virus invading.
- Thumb loop wrists prevent sleeves from sliding up when putting gloves on
- Color: blue
- Size: 125\*200cm, 62g
- Package: 1pc/bag, 100pcs/carton
- Carton size: 460\*200\*250mm

## YM-G016 LAB COAT

### Specification:

- Materials: Non-Woven PP/PE
- Color: White/blue/Yellow/Green
- Weight: 30/35/40/45GSM
- Size: M: 105\*130cm /L: 105\*135cm/XL: 115\*140cm/XXL: 115\*145cm
- With Label Collar /Knitted Collarand 4 Press Buttons Button Closure or Zipper Closure or Velcro Closure
- Size: M/L/XL/XXL

- Suited To Hospital and Daily Home Use Form of Packing: 10Pcs/Inner poly bag 100boxes/Carton
- Size of packing:
- Export carton: 560\*320\*420mm
- Gross Weight: 12.5kgs Net Weight: 11.0kgs



## YM-G017 LAB COAT

### Specification:

- Non Woven PP Fabric 40g, 45g
- 71 \*62cm Short sleeve Top
- 4Button Front Closure
- 108\*56cm Pant With Elastic Band
- Color: blue
- Keeping The Environment Free of infection and
- Disease is Always a Battle
- For a Health Professional

- Suited To Hospital and Daily Home Use
- Form of Packing: 1pcs/Inner poly bag
- 50boxes/Carton
- Size of packing:
- Export carton: 560\*350\*500mm
- Gross Weight: 12.5kgs Net Weight: 11.0kgs



## YM-G018 DISPOSABLE SCRUB COAT

### Specification:

- Disposable Scrub Coat are Perfect For any Medical Professional Shopping on a Budget. Available In Several Colors And Made From Quality Non Woven PP Or SMS Material, Disposable Scrub Coat Are Attractive as Well as Affordable Suited To Hospital And Daily Home Use
- Elastic Waist Pant
- Color: Green, blue
- Unisex
- Form of Packing: lpcs/Inner poly bag
- SOboxes/Car ton

### Description:

- Non Woven PP or SMS Fabric: 48g, 55g
- V Collar Top
- Size of packing:
- Export carton: 430\*340\*550mm
- Gross Weight: 8.5kgs Net Weight: 7.5kgs

## YM-G019 OVERALL GOWN

### Specification:

- With Head Hood and Shoe Cover and Zipper
- Materials: PP/PP+PE Film/SMS
- Color: White/blue/ Green/Dark
- Weight: 35/45/50/GOGSM,
- Size: M/L/XL/XXL at Clients'option
- S: 164\*126cm/ M: 172\*130cm /L: 178\*134cm/ XL: 182\*138cm 1
- XXL: 184\*142cm
- Suited To Lab, Industrial, Agriculture, Hospital And Daily Home Use
- Unique One Piece Design That Cover Whole Body Convenience And Daily Home Use.
- Unique One Piece Design That Cover Whole Body Convenience And
- Easy To Wear And Off
- Available in Standard or PE Coated
- Coverall, Zipper With Storm Flap Attached Hood,
- Boots Elastic Wrists, Waists and Ankles Are
- Available Bound Seam Construction
- Packing: lpcsjpoly bag
- SObags/Carton
- Carton size: 480\*340\*410mm





## YM-G020 STERILE DELIVERY PACK


- Non-woven surgical gown -IPC
- 86\*80cm Packing Clothing 2ply
- (Non-woven and thin film)-IPC
- 97\*73cm Sheet 3ply
- (thin film, Non-woven and paper pulp)-IPC
- 97\*73cm Triangle Trouser 3ply
- (thin film, Non-woven and paper pulp)-2PCS
- 66\*43cm Non-woven Cloth 2ply
- (thin film and Non-woven)-IPC
- 66\*43cm Non-woven Cloth 1ply-IPC;
- Umbilical Cord Clamp-IPC
- 5cm Cotton sticker-2PCS; PE glove-IPAIRS;
- Latex examination glove-IPAIRS
- 8\*10cm Gauze swab-5PCS; Face mask-IPC; Cap-IPC

## YM-G021 PE SHOE COVER

### Specification:

- Color: blue and Green
- Materials: PE Fabric: 1.6g, 1.8g, 2.0g, 2.1g, 2.5g,
- Made of PE With Comfortable Elastic Opening.
- Available in Different Thickness and color

- Size: 36\*15cm, 41 \*15cm, 40\*16cm
- Thickness: 0.0216mm, 0.018mm, 0.02mm, 0.021mm, 0.025mm at
- Clients'option
- Packing: 100pcs/poly bag
- 5000pcs/Carton
- Carton size: 480\*320\*360mm



## YM-G022 CPE SHOE COVER

### Specification:

- Color: blue and Green
- Materials: CPE Fabric: 1.8g, 2.0g, 2.1g, 2.5g, 3.0g
- Made of CPE Embossed Surface
- With Comfortable Elastic Opening
- Available in Different Thickness and color

- Size: 36\*15cm, 41 \*15cm, 40\* 16cm
- Thickness: 0.02mm, 0.025mm, 0.028mm, 0.035mm, at Clients'option
- Packing: 100pcs/poly bag
- 5000pcs/Carton
- Carton size: 480\*320\*360mm



## YM-G023 NON-WOVEN SHOE COVER

### Specification:

- Color: White and blue
- Materials: PP Non-woven
- Type A: PP Regular Shoe Covers
- Comfortable Elastic Hem, Helps You Keep Your Cool.
- Easily Worn Under A Sterile Gown. One-size-fits-most Sizing Reduces
- Your Inventory
- Type B: PP Non-skid Shoe Covers
- Non-skid Shoe Covers are Designed for Those Situations Where
- Wet or Slippery
- Conditions Exist
- A Proprietary Plastic Resin is Applied to The Bottom of The Cover
- to improve The Wearers Traction On The Floor
- Weight: 30/40/45gsm at Clients'option
- Size: 36\*15cm, 38\*18cm, 40\*16cm,
- Packing: 100pcs/bag
- 1000pcs/Carton
- From of packing:
- 100pcs/bag 1000pcs/Carton
- Size of packing: Export carton: 520\*280\*350mm

## YM-G024 PE SLEEVE

### Specification:

- Size: 20\*40cm, 42\*22cm, 46\*22cm
- Thickness 0.02 mm/0.018mm
- Made of PE .With Comfortable Elastic Water Repel len
- Elasticized at Both Ends to Prevent Slippage Fluid Resistant
- Effective Protection For Wrist and Forearm Clear, White, Blue, Red Color etc.
- Packing: 100pcs/poly bag
- 2000pcs/Carton
- Carton size: 340\*230\*320mm



## YM-G025 PE APRON

### Specification:

- Made of Polyethylene
- With Adjustable Tie
- Available In Different Size and Color
- Thickness 0.018 mm, 0.020mm, 0.025mm
- Size: 110\*70cm, 135\*71cm, 160\*90cm, etc.
- Color: White and blue
- Package: 100pcs/bag, 1000pcs/Carton
- Packing: 100pcs/poly bag, 1000pcs/Carton
- Carton size: 360\*250\*180mm




## YM-G026 UNDERPAD

- Material: Non-Woven fabrics+ Fluff pulp+SAP+PE Film
- Color: blue and green
- 60\*60cm 10pcs/bag 10bags/ctn 57\*22\*31cm
- 60\*90cm 10pcs/bag 10bags/ctn 42\*28\*42cm

## YM-G027 NON-WOVEN FABRICS BEDSHEET

### Description:

- Soft touch and harmless to your skin
- 100% biodegradable and environmentally friendly
- SPEC: Light blue, 60\*180cm
- 30gsm other size and specification are also acceptable
- PACKING: 25pcs/polybag, 4bags/ctn
- CTN SIZE: 32\*25\*24cm/100pcs



## YM-G028 BED SHEET

- Material: Film+paper
- Size: 90x190cm, etc.
- Color: Blue, green, etc.
- Packing: 10pcs/polybag, 100pcs/ctn

## YM-G029 ADULT DIAPER

### Specification:

- 65\*80cm, with leg cuffs, waistline: 71-121cm
- 65\*80cm, without leg cuffs, waistline: 71-111cm
- 80\*96cm, with leg cuffs, waistline: 96-147cm
- 80\*96cm, without leg cuffs, waistline: 96-137cm

### Material:

- Inner layer: Polyolefin Non-Woven
- Absorbent core, dry and comfortable.
- Refastening Frontal tape is convenient to put on and removed
- PE film makes the surface dry and comfort
- Super soft non-woven surface
- Leak Guard prevent leakage efficiently
- Package: 10pcs/bag, 80pcs/carton
- Carton size: 680\*255\*480mm

